



ULTRASMOOTH *Aqua*®

Tiemann Tip
Hydrophilic Coated Catheter
with Water Sachet



NEW
Tiemann
Tip

ULTRASMOOTH Aqua[®]

Tiemann Tip Hydrophilic Coated Catheter with Water Sachet

For **Complimentary Samples** simply call
FREEPHONE
0800 652 0424

Product Information Label

Self-adhesive product information label.



Notch on Funnel

Shows the alignment of the Tiemann Tip while being inserted.

Integral Sterile Water Sachet

To activate the hydrophilic coating. Easy to use, ready to go.

Green Introducer Grip

Soft but firm **green introducer grip** to assist with a more precise and hygienic insertion.

Non-coated grip zone.

Tiemann Tip

Curved, angled upwards to negotiate the male urethra.

Ultra Smooth Polished Eyes

Engineered to an exceptionally **smooth finish** to avoid the risk of urethral trauma.

Safe and Reliable Hydrophilic Coating

Activated in 15-30 seconds, remains smooth and highly lubricated from insertion to removal.



Dual Purpose Self-adhesive Tab

For easy catheter access and attachment to a convenient surface.



UltraSmooth Aqua® Tiemann Tip Ordering Codes

Male Length (42cm) – Pack of 30

Size	Order Code
● 10Ch	USAT4210
○ 12Ch	USAT4212
● 14Ch	USAT4214
● 16Ch	USAT4216

All **UltraSmooth Aqua® Tiemann Tip** catheters are available on prescription from your GP and can be obtained from your usual supplier.

For **Complimentary Samples** simply call
FREephone 0800 652 0424

For more information on our **UltraSmooth Aqua® Tiemann Tip** catheters, please visit our website:
www.peakmedical.co.uk

We are social

