



ULTRASMOOTH
Meatal Dilator™

UltraSmooth Meatal Dilator

Building on the successful **UltraSmooth** brand of intermittent catheters Peak Medical introduced the **UltraSmooth Meatal Dilator**. The Meatal Dilator is 12cm long and is a single use hydrophilic coated catheter without drainage eyes, designed specifically for dilatation of meatal strictures. Available in various Charrier sizes to suit all needs.

The benefits to you...

- Durable, consistently stable and highly lubricous, the hydrophilic coating ensures ultra smooth insertion and removal.
- The durable hydrophilic coating resists 'drying out', unlike other hydrophilic coatings.
- The coating is activated in water within 15 seconds.
- Finger holes allow primary packaging to be opened easily.
- Integrated funnel for grip and no-touch insertion.
- An adhesive pad on the packaging allows it to stick securely to most surfaces for preparation.
- Foil packaging ensures that the product reaches the end user in optimum condition.

Available on NHS prescription.

UltraSmooth Meatal Dilator Product Codes

Meatal Dilator Length 12cm – Pack of 25

Charriere Size	Product Code
08Ch	USMD08
10Ch	USMD10
12Ch	USMD12
14Ch	USMD14
16Ch	USMD16
18Ch	USMD18



For more information please call Freephone 0800 652 0424

