

Ileomate 1-Piece Drainable Pouches

The **Ileomate 1-Piece Drainable Pouch** has been designed to enhance comfort, whilst ensuring maximum security with the unique **SoftSafe** closure. Available in various capacities and sizes with the choice of cut to fit or pre-cut skin barriers.

The benefits to you...

Easier Drainage

- Unique **Finger Pocket** for easy drainage and hygienic cleaning of the pouch.
- Simply insert one or two fingers into the pocket.



Secure Closure

- Unique **SoftSafe** closure.
- Easy closure system.
- Providing optimum comfort with maximum security.



Skin Barrier

- Traditionally formulated, slim, flexible, skin friendly, hydrocolloid adhesive skin barrier that effectively absorbs moisture, maintaining the skin in optimum condition.
- Available in various shapes that offer maximum security in areas where creases and skin folds are likely to occur around the stoma.

Effective Bar Filter

- High capacity 32mm bar filter.
- Helps prevent ballooning and pancaking, contributing to improved pouch wear time.

Filter Protection

- Enhanced by **Gore –Tex®** membrane to prevent 'strike through' from liquid output.
- An additional layer of film, the 'raincoat', separates the filter from the stoma, providing optimum performance.



Pouch Cover

- Body side covering is made of a soft, non-woven material that helps to prevent skin soreness by absorbing moisture from the skin.
- Available in beige for discretion or in a 'needle pin' split cover to assist with pouch application and observation of the stoma.

This product range is latex free.

Available on NHS prescription.

Ileomate 1-Piece Drainable Pouches Product Codes

Small Drainable Pouch – Pack of 30

Stoma Size	Transparent	Beige
Cut to fit 13-60mm	ISTV113	ISBV113



Medium Drainable Pouch – Pack of 30

Stoma Size	Transparent	Beige
Cut to fit 13-60mm	IMTV113	IMBV113
Pre-cut 19mm	IMTV119	IMBV119
Pre-cut 22mm	IMTV122	IMBV122
Pre-cut 25mm	IMTV125	IMBV125
Pre-cut 28mm	IMTV128	IMBV128
Pre-cut 32mm	IMTV132	IMBV132



Large Drainable Pouch – Pack of 30

Stoma Size	Transparent	Beige
Cut to fit 13-60mm	ILTV113	ILBV113
Pre-cut 19mm	ILTV119	ILBV119
Pre-cut 22mm	ILTV122	ILBV122
Pre-cut 25mm	ILTV125	ILBV125
Pre-cut 28mm	ILTV128	ILBV128
Pre-cut 32mm	ILTV132	ILBV132
Pre-cut 35mm	ILTV135	ILBV135
Pre-cut 38mm	ILTV138	ILBV138
Pre-cut 45mm	ILTV145	ILBV145



Large Drainable Pouch with Oval Barrier – Pack of 30

Stoma Size	Transparent	Beige
Cut to fit 10-110mm	ILTOV110	-



Extra Large Drainable Pouch with Oval Barrier – Pack of 30

Stoma Size	Transparent	Beige
Cut to fit 10-110mm	IXLTV110	-

