




ULTRASMOOTH MCath4™

Non-coated Nelaton Catheter
with 4 Drainage Eyes
and Introducer Grip



ULTRASMMOOTH MCath4™

For Complimentary Samples simply call
FREephone 0800 652 0424

Non-coated Nelaton Catheter with
4 Drainage Eyes and Introducer Grip

Designed specifically for the drainage of urine from a
Mitrofanoff catheterisable channel. MCath4™ can also
be used for urethral catheterisation and ACE procedures.

**High Quality Smooth
Non-coated PVC Catheter**

Requiring water soluble
lubrication prior to use.

Flexible Funnel

For a secure connection to a
night drainage or other urine
collection system.

Ultra Smooth Polished Eyes


Four staggered heat polished drainage eyes for improved
drainage, especially where mucus is a problem.
Engineered to an exceptionally **smooth finish**
for a safer insertion and removal.

Green Introducer Grip

Soft but firm **green introducer grip** to assist
with a more precise and hygienic insertion.

Product Information Label

Self-adhesive product information
label for convenient patient record.



UltraSmooth MCath4™ Ordering Codes

Length (40cm) – Pack of 60

Size	Order Code
● 08Ch	USM4008
● 10Ch	USM4010
○ 12Ch	USM4012
● 14Ch	USM4014
● 16Ch	USM4016
● 18Ch	USM4018



All **UltraSmooth MCath4™** catheters are available on prescription from your GP and can be obtained from your usual supplier.

For **Complimentary Samples** simply call

FREephone 0800 652 0424

For more information on our **UltraSmooth MCath4™ catheters**, please visit our website:

www.peakmedical.co.uk

We are social

