



Fleximate

Fleximate 2-Piece Flanges

The **Fleximate 2-Piece Flanges** provide an alternative to, and are interchangeable with **Coloplast's Assura*** or **SenSura*** 2-piece ostomy system, providing users with more choice.

The benefits to you...

- The Fleximate 2-piece flanges are designed with softness, flexibility and comfort in mind. This has been achieved by laminating the hydrocolloid skin barrier with a flexible thin layer of polymeric foam.
- To provide comfort and protection for the stoma, the flanges have no plastic film which can potentially damage the stoma.
- For mild and severe stomal retraction, we offer a range of **shallow** and **deep** convexity flanges**.
- All Fleximate flanges will securely fit **Coloplast's Assura*** or **SenSura*** **40mm, 50mm** and **60mm** 2-piece pouches.
- Available in corresponding flange gasket sizes: **40mm, 50mm** and **60mm**. Also available in a wide selection of pre-cut sizes which can easily be stretched to the exact stoma size.

*Coloplast Assura and SenSura are registered trademarks of Coloplast A/S Denmark.

**All convex products should be used upon the advice and guidance of a Stoma Care Nurse.

These products are latex free.

Available on NHS prescription.

For more information please call Freephone 0800 652 0424



Fleximate 2-Piece Flanges Product Codes

SuperFlex Flange – Pack of 5

Stoma Size	40mm Flange	50mm Flange	60mm Flange
13mm starter hole	FFB4013	FFB5013	FFB6013
Pre-cut 22mm	FFB4022	-	-
Pre-cut 25mm	FFB4025	FFB5025	-
Pre-cut 28mm	FFB4028	FFB5028	-
Pre-cut 32mm	-	FFB5032	FFB6032
Pre-cut 35mm	-	FFB5035	FFB6035
Pre-cut 38mm	-	FFB5038	FFB6038
Pre-cut 41mm	-	-	FFB6041
Pre-cut 44mm	-	-	FFB6044
Pre-cut 47mm	-	-	FFB6047



SuperFlex Flange – Pack of 5

Stoma Size	40mm Flange	50mm Flange	60mm Flange
13mm starter hole	FSS4013	FSS5013	FSS6013
Pre-cut 22mm	FSS4022	-	-
Pre-cut 25mm	FSS4025	FSS5025	-
Pre-cut 28mm	FSS4028	FSS5028	-
Pre-cut 32mm	-	FSS5032	FSS6032
Pre-cut 35mm	-	FSS5035	FSS6035
Pre-cut 38mm	-	FSS5038	FSS6038
Pre-cut 41mm	-	-	FSS6041
Pre-cut 44mm	-	-	FSS6044
Pre-cut 47mm	-	-	FSS6047



SuperFlex Flange With Shallow Convexity – Pack of 5

Stoma Size	40mm Flange	50mm Flange	60mm Flange
13mm starter hole	FC4013	FC5013	FC6013
Pre-cut 16mm	FC4016	-	-
Pre-cut 19mm	FC4019	-	-
Pre-cut 22mm	-	FC5022	-
Pre-cut 25mm	-	FC5025	-
Pre-cut 28mm	-	FC5028	-
Pre-cut 32mm	-	-	FC6032
Pre-cut 35mm	-	-	FC6035
Pre-cut 38mm	-	-	FC6038
Pre-cut 41mm	-	-	FC6041



SuperFlex Flange With Deep Convexity – Pack of 5

Stoma Size	50mm Flange	60mm Flange
13mm starter hole	FXC5013	FXC6013
Pre-cut 19mm	FXC5019	-
Pre-cut 22mm	FXC5022	-
Pre-cut 25mm	-	FXC6025
Pre-cut 28mm	-	FXC6028
Pre-cut 32mm	-	FXC6032



NEW

Fleximate Soft Convex Flanges Product Codes

Soft Convex Flange – Pack of 5

Stoma Size

13mm start hole
Pre-cut 16mm
Pre-cut 19mm
Pre-cut 22mm
Pre-cut 25mm
Pre-cut 30mm
Pre-cut 35mm
Pre-cut 41mm

60mm Flange

FSC6013-36
FSC6016
FSC6019
FSC6022
FSC6025
FSC6030
FSC6035
FSC6041



For more information please call Freephone 0800 652 0424