



ULTRASMOOTH PLUS™

NEW

UltraSmooth Plus ISC Catheters

The **UltraSmooth Plus** hydrophilic coated nelaton catheter incorporates the latest technologies to produce ultra smooth heat polished and coated drainage eyes, providing optimum smoothness and comfort. Available in three different lengths: 40cm (male and/or Mitrofanoff), 30cm (paediatric) and 18cm (female). They are available in various Charrier sizes to suit all needs.

The benefits to you...

- A soft but firm **green introducer grip** is provided to assist with a more hygienic insertion and removal (available with the male and paediatric catheter only).
- Durable, consistently stable and highly lubricious, the hydrophilic coating ensures ultra smooth insertion and removal.
- The durable hydrophilic coating resists 'drying out' unlike other hydrophilic coatings.
- The coating is activated in water within 15 seconds.
- 'Finger holes' allow primary packaging to be opened easily.
- Integrated funnel for grip and no-touch insertion. Allows connection to a variety of external collection devices.
- Adhesive pad on packaging allows it to stick securely to most surfaces for preparation.
- Self-adhesive product information label for convenient patient record.
- Foil packaging ensures that the product reaches the end user in optimum condition.

Available on NHS prescription.

For more information please call Freephone 0800 652 0424



UltraSmooth Plus ISC Catheters

Product Codes

Male Length 40cm – Pack of 25

Carrier Size	Product Code
08Ch	USP4008
10Ch	USP4010
12Ch	USP4012
14Ch	USP4014
16Ch	USP4016
18Ch	USP4018



Female Length 18cm – Pack of 25

Carrier Size	Product Code
08Ch	USP1808
10Ch	USP1810
12Ch	USP1812
14Ch	USP1814
16Ch	USP1816
18Ch	USP1818



Paediatric Length 30cm – Pack of 25

Carrier Size	Product Code
06Ch	USP3006
08Ch	USP3008
10Ch	USP3010




Self-Adhesive Sticky Pad

Assists preparation. Enables catheter packaging to be fixed to a convenient surface.



Product Information Label

Self-adhesive product information label for convenient patient record.



Easy Access Finger Holes

Require minimum effort when opening catheter packaging. Ideal for those with limited dexterity.