



Varimate Wedges

The innovative **Varimate Wedges** provide a simpler and safer way of managing problematic stomas with the user in mind.


The benefits to you...

- These double-sided hydrocolloid pre-cut wedges can be used for filling creases and folds around the stoma where leakage can occur.
- Easier to manage than conventional skin pastes and fillers.
- The wedges can easily be peeled off the 'wedge wheel' and stacked one on top of the other until the desired depth is achieved to create a platform for the appliance.
- Can also be used as an alternative to products with convexity.

Available on NHS prescription.

Varimate Wedges Product Codes

Size	Pack Size	Product Code
Small	60	VMW60
Large	40	VMW40



For more information please call Freephone 0800 652 0424

