



Uromate

## Uromate 1-Piece Urostomy Pouches

The **Uromate 1-Piece Urostomy Pouch** has been designed to enhance comfort, whilst ensuring maximum security. Available in various capacities and sizes with the choice of cut to fit or pre-cut skin barriers.

### The benefits to you...

#### Skin Barrier

- Traditionally formulated, slim, flexible, skin friendly, hydrocolloid adhesive skin barrier that effectively absorbs moisture, maintaining the skin in optimum condition.
- Available in various shapes that offer maximum security in areas where creases and skin folds are likely to occur around the stoma.

#### Pouch Cover

- Body side covering is made of a soft, non-woven material that helps to prevent skin soreness by absorbing moisture from the skin.
- Available in a 'needle pin' split cover to assist with pouch application and observation of the stoma.

#### Non-Return Valve

- A non-return valve prevents backflow of urine onto the stoma and skin barrier, preventing skin soreness and increasing pouch wear time.

#### Easier Drainage

- The wide bore, soft silicone bottom outlet tap allows drainage of urine and mucus.
- Wide bore adaptors included to allow connection to universal leg bags and overnight drainage systems.
- Unique '**Tuck-Away**' pocket for retaining the outlet tap for enhanced comfort and discretion.

**This product range is latex free.**

**Available on NHS prescription.**

For more information please call Freephone 0800 652 0424



# Uromate 1-Piece Urostomy Pouches Product Codes

## Large Urostomy Pouch – Pack of 20

| Stoma Size         | Transparent |
|--------------------|-------------|
| Cut to fit 13-60mm | ULT113B     |
| Pre-cut 19mm       | ULT119B     |
| Pre-cut 22mm       | ULT122B     |
| Pre-cut 25mm       | ULT125B     |
| Pre-cut 28mm       | ULT128B     |
| Pre-cut 32mm       | ULT132B     |
| Pre-cut 35mm       | ULT135B     |
| Pre-cut 38mm       | ULT138B     |



## Medium Urostomy Pouch – Pack of 20

| Stoma Size         | Transparent |
|--------------------|-------------|
| Cut to fit 13-60mm | UMT113B     |



## Extra Small Urostomy Pouch – Pack of 20

| Stoma Size         | Transparent |
|--------------------|-------------|
| Cut to fit 10-50mm | UXST110B    |



## Urostomy Night Drainage Adaptors – Pack of 10

**Product Code**  
AD10



## Night Drainage Bag – Pack of 10

**Product Code**  
NDB10

